# OPED

**OPED GmbH**  
Medizinpark 1  
83626 Valley/Oberlaindern  
Germany  
Fon +49 (0) 80 24/60 81 8-210  
Fax +49 (0) 80 24/60 81 8-299  
mail@oped.de  
www.oped.de

**OPED AG**  
Hinterbergstrasse 26  
6330 Cham  
Switzerland  
Fon +41 (0) 41/748 53 35  
Fax +41 (0) 41/748 53 40  
contact@oped.ch  
www.oped.ch

**OPED UK Ltd**  
Unit 5  
Avro Business Centre  
Avro Way  
Bowerhill  
Melksham  
SN12 6TP  
United Kingdom  
Phone +44 (0) 1380 722177  
Fax +44 (0) 1380 710095  
enquiries@oped.biz  
www.opeduk.co.uk

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Ihr OPED Team

## CAMOped

Intended purpose: CAMOped is a treatment device and is used to move the knee or the hip gently at an early stage of treatment.

### Safety Advice

- ⚠ The patient may only make adjustments in consultation with his/her treating doctor.
- ⚠ Should you experience excessive pain, swelling, numbness, sensitivity changes, unusual reactions or an uncomfortable sensation during use, please contact your doctor or orthopaedic technician immediately.
- ⚠ This treatment device was developed for the indications stated and may be damaged through unusual and non-intended use.
- ⚠ Please consult your treating doctor before using it in conjunction with another medical device.
- ⚠ The product must be cleaned and disinfected and its functions must be tested after every use.
- ⚠ The treatment device must undergo maintenance (disassembly, cleaning and functional tests) at least once per year.
- ⚠ Only replacement parts approved by OPEd GmbH may be used.
- ⚠ Note: Not made of natural rubber latex. 
- ⚠ In the event of malfunctions, repair work must be carried out by a technician trained and certified by OPEd GmbH.

### Assembling and dismantling the footrests

#### Assembling the footrests

1. Insert the right footrest into the right slider (marked in green). The black bolt must be placed between the orange adjustment caps. (a)
2. Insert the left footrest into the left slider (marked in red). The black bolt must be placed between the orange adjustment caps. (b)
3. Please note that the black knob must be open if you want to be able to move the footrest. (f)
4. Press the footrests down firmly with both hands until you hear a 'click'.

#### Marking on the device

- The right footrest and slider are marked in green.
- The left footrest and slider are marked in red.

#### Dismantling the footrests

1. Press down the small round button on the slider and remove the footrest. (c)
2. First place the CAMOped and then the boxes with the footrests into the cardboard box. (e)
3. When putting them back in the packaging, please ensure that both sliders are placed in the middle. Otherwise, the inserts for the frames won't fit perfectly.

### Preparation

#### Angle restriction

Use the angle restriction to set the permitted range that your leg can bend or stretch. Your treating doctor will have told you how far you are allowed to bend or stretch. If you are unsure, please consult your doctor/physiotherapist.

- Always set the left and right sliders the same on both sides (d)
- The numbers 0-12 roughly correspond to angles of 0°-120°

### How to Use – Five Steps to Success

#### 1. Position product (Fig 1)

- Place the CAMOped with extended supports between your legs.
- There should be a gap of approx. 10cm (one hand's width) between the device and your buttocks.

Tip: An instruction video for using on bed/couch/floor can be found at [www.camoped.com](http://www.camoped.com).

#### 2. Insert legs (Fig 2)

- Place your legs into the frames. Be careful when inserting the injured leg!
- Both legs should be slightly bent here and the sliders should be at the same height.
- The entire sole of the foot should have contact with the foot plate.

#### 3. Lower leg length adjustment (Fig 3)

- Adjust the frames individually by sliding them.

#### 4. Close the belt straps (Fig 4)

- Close all of the belt straps.

#### 5. Position the sliders (Fig 5)

Check whether the angle restrictions have been correctly set and snapped closed. (d)

- Press down the long, curved button and keep it held down to unlock the sliders.
- Extend one leg to the maximum extension permitted and bend the other leg to the maximum flexion angle allowed.
- Release the button again and the sliders will snap into the desired position.

#### 6. Start movement exercise (Fig 6)

- Start slowly and concentrate on the movement exercises.

## Product Characteristics

CAMOPed is available in one standard size.

## Indications

- Cruciate ligament ruptures
- Corrective osteotomies
- Knee and hip prostheses
- Reconstructive operations on the menisci and joint cartilage
- Following hip/knee arthroscopy in combination with synovectomy
- Retropatellar damage

## Contraindications

None known

## Technical data/parameters

### ● Materials used

Profile: Aluminium

Metal parts: Aluminium, steel

Plastic parts: PP, PA, POM, TPE, ABS

Frame padding: PE foam

### ● Dimensions (assembled)

83 cm x 60 cm x 40 cm

### ● Weight: 8.780 g

## Transport conditions/storage conditions/disposal:

**Storage:** Store in a dry place and out of direct sunlight.

## Guarantee

The guarantee is in line with statutory prescriptions.

## Cleaning and Care

- Strap: Machine wash at 40°C.
- Clean plastic parts and pads with a damp cloth and with plastic-friendly cleaning agents only.
- Storage: store in a dry place.

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